

GAS CABINET 3 CYL



- Sikkerhet gass kabinett
- Bruk med tre L50
- Egnet for rene gasser og giftige blandinger
- Laget av malt stål

Laboratorier og analyse

Domains of applications

Designed for implementation of:

Rene gasser og blandinger, giftige blandinger


Dimensions

Length (L) : 1100 mm

Height (H) : 2000 mm

Depth (D) : 525 mm

Net weight : 155 kg



Products

| Product Number | Short Designation |
|----------------|-------------------------|
| 8063 | GAS CABINET 3 CYLINDERS |