

RACK IND WITH CHAIN



- Gass-sylinder festesystem
- vegg~~POS=TRUNC montert~~POS=HEADCOMP
- Laget av malt stål

sveising, cutting, Varmebehandling, Industrial Process

Domains of applications

Designed for:


L50 og M20 komprimert gassflaske, B84 flytende gass sylind-
deren

Dimensions

Length (L) : 450 mm

Height (H) : 243 mm

Net weight : 1.27 kg



Products

Product Number	Short Designation	Material
162192	CYLINDER RACK 1.CYL UNIVERSAL WITH CHAIN	Karbonstål